



The name you can trust

Walk-In Humidity (Stability) Chambers



Green Mode Energy Optimisation



Low Power Consumption



More efficiency



Less breakdown



Less deviations



More storage area





ECOtronic Walk-in Stability Chambers for a greener tomorrow

Ecotronic Working Principle

Dewpoint Control

- ◆ Measure of atmospheric moisture
- ◆ Temperature to which air must be cooled in order to be saturated with water vapor already in the air
- ◆ Temperature at which dew or frost will form
- ◆ Coil temperature controlled as per Dewpoint and Psychometric chart



Features

- ◆ Works on the principle of Dewpoint Control
- ◆ No deviations during short-term door opening & power failure
- ◆ No utility Area | More storage area
- ◆ Lower operating cost | Minimum power & water consumption
- ◆ Less deviations | Less breakdowns
- ◆ Increase in system & process efficiencies through reduced energy consumption



Our Products Range: Humidity (Stability) Chambers | Photostability Chambers | Cooling Chambers | Walk-in Chambers | Deep Freezers | Plasma Freezers(-40°C) | Ultra Low Temperature Freezers(-80°C) | Incubators



The name you can trust

Visit us on: www.newtronic.in

NEWTRONIC LIFECARE EQUIPMENT PVT. LTD.

H.O.: 108-ABCD, Kandivali Co-op. Ind. Estate Ltd.
Charkop, Kandivali(West), Mumbai - 400067, India
T: +91 22 28679326 | F: +91 22 28676059

Email:

India : sales@newtronic.in

International : exports@newtronic.in

Visit us on: www.newtronic.in