



**Controlling humidity is always
great challenge !!!**

**We understand stability at its core...
One of the largest manufacturer & suppliers
of stability chambers with wide selection range**



200 ltr | 7 Cu Feet




430 ltr | 16 Cu Feet



730 ltr | 54 Cu Feet



1000 ltr



1500-2000 ltr. (Double Door)

Humidity (Stability) Chambers


Temperature Range: 20°C to 60°C
Humidity Range: 40% Rh to 95% Rh

We offer a centralized solution across all ICH conditions of real-time, intermediate, accelerated & semi-permeable study conditions along with a 21 CFR Part 11 ready software that provides tools for viewing, trending, alarm management, audit trails, MKT & many other such features.


Salient features




Tech Specs



Telescopic Tray



High Quality Temp.-RH Sensor



Forced Air Circulation



Standby Humidity with booster humidifier System

Innovative features



Alarms



Safety devices



Remote Alarm Contact



ICDAS Compatible



MEGVID Compatible


Humidity (Stability) Chambers

Temperature Range: 20°C to 60°C
Humidity Range: 40% Rh to 95% Rh

- We are one stop solutions for stability requirement that includes:
- Installations and Commissioning
- DQ, IQ, OQ, PQ documentation
- Training
- Online Support
- Annual Maintenance Contract


Applications

Newtronic Chambers are suitable to meet the regulatory requirements of the various International bodies such as ICH, FDA, TAPPI, TGA, JPMMA, MHRA, ANVISA, ASEAN and other testing standards




- General Conditions
- Semi-permeable Conditions
- MVTR and Other Special Conditions

Industrial Applications



Pharmaceutical



Biotechnology



Cosmetic



Food & Beverages

SPECIFICATIONS

- Temperature Range : 20°C to 60°C
- Humidity Range : 40% Rh to 95% Rh
- Temperature Accuracy : ± 0.2°C
- Humidity Accuracy : ± 2% Rh
- Temperature Uniformity : ± 1°C
- Humidity Uniformity : ± 3% Rh

Power supply 230 Volts 50 Hz main supply. Also available with 60 Hz cycle

Optional Accessories

- ICDAS - A Single Software installation for all Newtronic make equipments
- 4 pt. Temperature & Humidity Scanner
- Stand by Cooling System with auto switch over (430 litre & above)
- Door Access system (Numeric or Biometric)
- Mobile Alert system
- Remote alarm system suitable for upto 4 nos. Newtronic make equipments

Selection Guide

Model	Capacity		Internal Size H* w* D*(mm)	External Size H* w* D*(mm)	Number of Trays	Tray Size (mm)	Equipment Weight (kg)
	Cu Ft.	Liter					
NLHC7SI	7	200	800 x 500 x 500	1295 x 640 x 990	3	430 X 460	150
NLHC16SI	16	430	1200 x 600 x 600	1695 x 740 x 1110	4	530 X 560	250
NLHC26SI	26	735	1500 x 700 x 700	2050 x 840 x 1210	5	630 X 660	275
NLHC34SI	34	1000	1500 x 800 x 800	2050 x 940 x 1310	5	730 X 760	350
NLHC54SI	54	1500	1500 x 1600 x 600	2050 x 1740 x 1125	10	695 X 560	350
NLHC71SI	71	2000	1500 x 1600 x 800	2050 x 1740 x 1325	10	695 X 760	375

* Newtronic reserve the rights to change the product specifications without prior notice